NADAL® PSA Test (test cassette)

REF 602003

DE	Gebrauchsanweisung	2	CZ	Návod k použití	30
EN	Instructions for use	6	FI	Käyttöohje	34
FR	Instructions d'utilisation	10		Symbols	39
ES	Instrucciones de uso	14		Our Teams	40
IT	Istruzioni per l'uso	18			
PL	Sposób użycia	22			
PT	Instruções de Utilização	26			

IVD **C** **€0197**

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2. Ballentine Carter H. (2004) Prostate Cancers in Men with Low PSA Levels – Must We Find Them? N Engl J Med 350(22): 2292-2294
3. Jain S et al. (2002) Improving the utility of prostate specific antigen (PSA) in the diagnosis of prostate cancer: the use of PSA derivatives and novel markers. Postgrad Med J (78) 646-650
4. Price et al. (2001) Pre- and post-analytical factors that may influence use of serum prostate specific antigen and its isoforms in a screening programme for prostate cancer. Ann Clin Biochem (38) 188-216
5. Ferlay J et al. (2013) Cancer incidence and mortality patterns in Europe: Estimates for 40 countries in 2012. Europ J Cancer (49) 1374-1403
6. Leitlinienprogramm Onkologie (Deutsche Krebsgesellschaft, Deutsche Krebshilfe, AWMF): Interdisziplinäre Leitlinie der Qualität S3 zur Früherkennung, Diagnose und Therapie der verschiedenen Stadien des Prostatakarzinoms, Langversion 3.0, 2014 AWMF Registernummer: 034/022OL, <http://leitlinienprogramm-onkologie.de/Leitlinien.7.0.html> (Zugriff am: 19.11.2015)

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Rev. 1, 2016-05-10 OM/JBo

1. Intended Use

The NADAL® PSA Test is a rapid visual immunoassay for the semiquantitative presumptive detection of prostate-specific antigen (PSA) in human serum, plasma or whole blood specimens with a cut-off of 4 ng/ml. The test detects total PSA (tPSA). The NADAL® PSA Test is intended for use as an aid in the diagnosis of prostate cancer by professional users as elevated PSA levels frequently indicate an increased risk of prostate carcinomas.

2. Introduction and Clinical Significances

Prostate cancer is the most frequent type of cancer found in men and the third most common cause of death due to cancer in men in Germany. Age-standardised incidences of prostate cancer differ worldwide. In Europe there is a decline from North to South. In the USA, incidences of prostate cancer are the highest, with Afro-American males especially affected. Incidences increase dramatically in males over 40 years of age, occurring in 50% of those over 70 years. The introduction of prostate specific antigen (PSA) as a tumor marker has led to a significant increase in the detection of prostate carcinomas, especially those in early stages. Looking at PSA from a biological viewpoint, it is a 33 kDa protein that is synthesised in the prostatic gland. It functions as a serine protease and serves to liquefy the seminal fluid. It is present in normal, benign hyperplastic and malignant prostate tissues as well as in metastatic prostate cancer and in seminal fluid. An elevation of the serum concentration is reported in patients with both benign prostatic hypertrophy (BPH) and prostate carcinoma, but rarely in healthy men. PSA is almost absent in the serum of healthy women. The PSA level in serum or plasma of normal healthy men should be lower than 4 ng/ml. If the structural integrity of the prostate is disturbed and/or the gland size is increased, the amount of PSA in blood/serum/plasma may rise, reaching levels of up to 200 ng/ml. At the cut-off of 4 ng/ml PSA, further medical analysis is recommended, even though in the PSA concentration range of 4-10 ng/ml the elevated levels are frequently not caused by cancer but by other factors, such as benign prostatic hyperplasia or prostatitis. PSA concentrations of >10 ng/ml strongly indicate the presence of prostatic carcinoma. Although an ethnical and/or age-specific modification of the cut-off has been discussed in the literature, the PSA amount of 4 ng/ml is the generally accepted value at which follow-up examinations of the patient should be carried out.

3. Test Principle

The NADAL® PSA Test enables the detection of prostate-specific antigen through visual interpretation of colour development on the internal test strip. Anti-PSA antibodies are immobilised in the test line region of the membrane. During testing, the specimen reacts with further anti-PSA antibodies conjugated to coloured particles and precoated onto the sample pad of the internal test strip. The mixture then migrates along the membrane by capillary action and interacts with reagents on the membrane. If the PSA concentration in the sample is ≥ 4 ng/ml, a coloured line appears in the test line region (T) of the membrane. The reference line (R) forms independently of the presence of PSA. At PSA concentrations of around 10 ng/ml the colour intensity

of the test line (T) and that of the reference line (R) match allowing a semiquantitative estimation of the PSA level. The appearance of a coloured line in the control line region (C) serves as a procedural control, indicating that a proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 20 NADAL® PSA test cassettes (incl. disposable pipettes)
- 1 Buffer
- 1 package insert

5. Additional Materials Required

- Timer
- Specimen collection container
- Centrifuge (for serum/plasma only)
- Lancets (for fingerstick whole blood only)

6. Storage & Stability

The test kit should be stored at 2-30°C until the expiry date printed on the packaging. The test cassette must remain in the sealed foil pouch until use. Do not freeze tests. Care should be taken to protect the components of the test kit from contamination. Do not use tests if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipments, containers or reagents can lead to false results.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure prior to testing.
- Do not use the test beyond the expiration date indicated on the package.
- Do not use the test if the foil pouch is damaged.
- Do not reuse tests.
- Do not add samples to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).
- Avoid cross-contamination of specimens by using a new specimen collection container for each specimen obtained.
- Do not substitute or mix components from different test kits.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being assayed.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious, and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Humidity and temperature can adversely affect test results.
- Used testing materials should be discarded according to local regulations.

8. Specimen Collection and Preparation

- The NADAL® PSA Test is intended for use with human serum, plasma or whole blood specimens only.
- Intake of Finasteride (5-alpha-reductase inhibitor) reduces the PSA concentration by a maximum of 50%. This should be considered in the interpretation of the results.
- Different factors may increase the PSA level and should be avoided during the following intervals before sample collection. Please instruct your patients in this regard:
 - Cycling: 24 hours before sample collection.
 - Sexual activities (ejaculation): 24-48 hours before sample collection.
 - Manipulation of the prostate in medical examinations. The following intervals are recommended before sample collection:

Examination	Interval
Prostate biopsy	> 6 weeks
Transurethral resection of the prostate	> 6 weeks
Transrectal prostatic ultrasound	> 1 week
Rigid cystoscopy	> 1 week
Digital rectal examination	3 days – 1 week
Prostate massage	> 1 week

- Containers containing anticoagulants such as EDTA, citrate or heparin should be used for plasma preparation or venous whole blood samples. Only clear, non-hemolysed specimens are recommended for use with this test. Serum or plasma should be separated from blood as soon as possible to avoid hemolysis.
- Perform testing immediately after specimen collection. Do not leave specimens at room temperature for prolonged periods. Serum and plasma specimens can be stored at 2-8°C for up to 3 days. For long-term storage, specimens should be kept below -20°C. Whole blood collected by venipuncture should be stored at 2-8°C if the test is to be run within 2 days of collection. Do not freeze whole blood specimens. Whole blood collected by fingerstick should be tested immediately.
- Bring specimens to room temperature prior to testing. Frozen specimens should be completely thawed and mixed well prior to testing. Avoid repeated freezing and thawing of specimens.
- If specimens are to be shipped, they should be packed in compliance with all applicable regulations for transportation of etiological agents.
- Icteric, lipemic, hemolysed, heat treated and contaminated specimens may cause erroneous results.

9. Test Procedure

Bring tests, specimens, buffer and/or controls to room temperature (15-30°C) before use. Do not open foil pouches until you are ready to perform the assay.

1. Remove the test cassette from the foil pouch and place it on a clean, level surface. Label the test cassette with a patient or control identification. The best results will be obtained, if the assay is performed within one hour after opening the foil pouch.

2. a) For serum or plasma specimens:

Holding the provided disposable pipette vertically, transfer 1 drop of serum/plasma to the specimen well (S) of the test cassette.

b) For venipuncture whole blood specimens:

Holding the provided disposable pipette vertically, transfer 2 drops of whole blood to the specimen well (S) of the test cassette.

c) For fingerstick whole blood specimens:

Allow 2 hanging drops of fingerstick whole blood to fall into the center of the specimen well (S) of the test cassette.

Avoid trapping air bubbles in the specimen well (S) and do not add any solution to the result area.

As the test begins to run, you will observe a coloured liquid migrate along the membrane.

3. Add 1 drop of buffer to the sample well (S). Start the timer.
4. Wait for the coloured lines to appear. The result should be read after 5 minutes. Do not interpret the result after more than 10 minutes.

10. Result Interpretation

Positive

Three coloured lines appear on the membrane. One line appears in the control line region (C) and the other two lines appear respectively in the test line region (T) and the reference line region (R).

- If the colour intensity of the test line (T) is weaker than that of the reference line (R), this indicates a PSA level of between 4 ng/ml and 10 ng/ml.
- If the colour intensity of the test line (T) is equal or close to that of the reference line (R), this indicates a PSA level of approximately 10 ng/ml.
- If the colour intensity of the test line (T) is stronger than that of the reference line (R), this indicates a PSA level of above 10 ng/ml.

Negative

Two coloured lines appear on the membrane. One line appears in the control line region (C) and the other line appears in the reference line region (R). No coloured line appears in

the test line region (T). This result indicates that the PSA concentration in the sample is less than 4ng/ml.

Invalid

The control line (C) and/or the reference line (R) fail(s) to appear. Results from any test which has not produced each of these lines at the specified reading time should be discarded. Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Note:

- At the cut-off concentration of 4 ng/ml, the colour intensity of the test line (T) is faint. Therefore, any shade of colour in the test line region (T) should be considered positive.
- Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line failure.

11. Quality Control

An internal procedural control is included in the test cassette:

A coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Good laboratory practice (GLP) recommends the use of control materials to ensure proper test kit performance.

12. Limitations

- The NADAL® PSA Test is for professional *in-vitro* diagnostic use and should only be used for the semiquantitative detection of PSA in human serum, plasma or whole blood specimens.
- The NADAL® PSA Test only indicates the presence of PSA in the specimen and should not be used as the sole criterion for the diagnosis of prostate cancer. The test is intended for use only as an aid in diagnosis of prostate cancer, as elevated PSA levels are only an indicator of an increased risk of prostate cancer.
- A significant number of patients with BPH (over 15%) and fewer than 1% of healthy individuals have elevated PSA levels. As with all diagnostic tests, a confirmed diagnosis should only be made by a physician after all clinical and laboratory findings have been evaluated.
- The cut-off of 4 ng/ml is widely accepted; however, age-specific and ethnic-specific modifications have been discussed. Before using the assay, ensure that 4 ng/ml is a suitable cut-off for your purpose, e.g. for men younger than 45 years this cut-off might be too high.
- It is recommended to check PSA levels at regular intervals, even in light of negative test results.
- The test should not be used for monitoring PSA levels after prostatectomy, as the expected levels are too low to be detected with the assay.
- Line intensities should not be used to quantify PSA concentrations. The test result is semiquantitative and only allows rough estimation of concentration ranges of

≥4 ng/ml and <10 ng/ml, ≈10 ng/ml, >10 ng/ml. This resolution is not sufficient to monitor the increase of PSA concentrations over the time.

- Negative test results should not be used to exclude the risk of prostate cancer as some tumors do not lead to elevated PSA levels. Therefore, the test is not meant to substitute conventional diagnostic examinations but rather to complement them. By this, PSA testing increases the detection of prostate carcinomas in their early stages.
- PSA concentrations can be transiently elevated or decreased by conditions listed under point 8 'Specimen Collection and Preparation'. These conditions should be taken into account when making the diagnosis.
- Additional individual risk assessments for prostate cancer (e.g. quantitative PSA testing, PSA density, PSA velocity or ratio of free to complexed PSA, prostate imaging) may be helpful to decide which follow-up measures should be taken for patients with positive test results.
- The test reliably detects 4 ng/ml over the appearance of a faint test line (T). In order to guarantee reliable test line formation at 4 ng/ml, it cannot be ruled out that concentrations slightly under 4 ng/ml may be occasionally detected with very faint lines.
- Specimens from patients who have received monoclonal antibodies from mouse for diagnostic or therapeutic purposes may contain human anti-mouse antibodies (HAMA). Such specimens may show either elevated or lowered values when tested with test kits that utilise monoclonal antibodies from mouse.

13. Performance Characteristics

Calibration

The test was calibrated against the international standard of the WHO (Prostate-Specific Antigen Free, NIBSC code: 96/668).

Clinical performance

244 positive samples (PSA ≥ 4 ng/ml) and 361 negative samples (PSA < 4 ng/ml) were determined with an EIA and tested with the NADAL® PSA Test. The relative diagnostic performance for the cut-off of 4ng/ml is summarised in the following table:

The NADAL® PSA Test vs. EIA

		EIA		Total
		+	-	
NADAL® PSA Test	+	241	5	246
	-	3	356	359
	Total	244	361	605

Relative sensitivity: 98.8% (97.3%-100.0%)

Relative specificity: 98.6% (97.6%-100.0%)

Overall agreement: 98.7% (98.7%-98.8%)

For the accuracy of the semiquantitative estimation of concentration ranges the following results were obtained:

	Concentration range [ng/ml]	EIA			Total
		<4	4-10	>10	
NADAL® PSA Test	<4	356	3	0	359
	4-10	5	86	3	94
	>10	0	4	148	152
	Total	361	93	151	605

Agreement (<4): 98.6%

Agreement (4-10): 92.5%

Agreement (>10): 98.0%

Interference testing

The following substances did not interfere with the test at the listed concentrations:

Ascorbic acid, gentisic acid, acetaminophen, acetylsalicylic acid, caffeine, uric acid (each 20 mg/dl), hemoglobin (1,000 mg/dl), oxalic acid (60 mg/dl), bilirubin (1,000 mg/dl), triglycerides (3,000 mg/dl).

Intra- and Inter-LOT variability

PSA concentrations of 0, 2, 4, 10, and 20 ng/ml were determined correctly > 99.9% when tested in 10-fold determinations with three independent LOTs.

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14. Bibliographie

- Milford Ward A. et al. (2001) Prostate specific antigen: biology, biochemistry and available commercial assays. Ann Clin Biochem 38: 633-651

- Ballentine Carter H. (2004) Prostate Cancers in Men with Low PSA Levels – Must We Find Them? N Engl J Med 350(22): 2292-2294
- Jain S et al. (2002) Improving the utility of prostate specific antigen (PSA) in the diagnosis of prostate cancer: the use of PSA derivatives and novel markers. Postgrad Med J (78) 646-650
- Price et al. (2001) Pre- and post-analytical factors that may influence use of serum prostate specific antigen and its isoforms in a screening programme for prostate cancer. Ann Clin Biochem (38) 188-216
- Ferlay J et al. (2013) Cancer incidence and mortality patterns in Europe: Estimates for 40 countries in 2012. Europ J Cancer (49) 1374-1403
- Leitlinienprogramm Onkologie (Deutsche Krebsgesellschaft, Deutsche Krebshilfe, AWMF): Interdisziplinäre Leitlinie der Qualität S3 zur Früherkennung, Diagnose und Therapie der verschiedenen Stadien des Prostatakarzinoms, Langversion 3.0, 2014 AWMF Registernummer: 034/022OL, <http://leitlinienprogramm-onkologie.de/Leitlinien.7.0.html> (Zugriff am: 19.11.2015)

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14. Referencias

1. Milford Ward A. et al. (2001) Prostate specific antigen: biology, biochemistry and available commercial assays. Ann Clin Biochem 38: 633-651
2. Ballentine Carter H. (2004) Prostate Cancers in Men with Low PSA Levels – Must We Find Them? N Engl J Med 350(22): 2292-2294
3. Jain S et al. (2002) Improving the utility of prostate specific antigen (PSA) in the diagnosis of prostate cancer: the use of PSA derivatives and novel markers. Postgrad Med J (78) 646-650
4. Price et al. (2001) Pre- and post-analytical factors that may influence use of serum prostate specific antigen and its isoforms in a screening programme for prostate cancer. Ann Clin Biochem (38) 188-216
5. Ferlay J et al. (2013) Cancer incidence and mortality patterns in Europe: Estimates for 40 countries in 2012. Europ J Cancer (49) 1374-1403

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14. Bibliografia

- Milford Ward A. et al. (2001) Prostate specific antigen: biology, biochemistry and available commercial assays. Ann Clin Biochem 38: 633-651
- Ballentine Carter H. (2004) Prostate Cancers in Men with Low PSA Levels – Must We Find Them? N Engl J Med 350(22): 2292-2294
- Jain S et al. (2002) Improving the utility of prostate specific antigen (PSA) in the diagnosis of prostate cancer: the use of PSA derivatives and novel markers. Postgrad Med J (78) 646-650
- Price et al. (2001) Pre- and post-analytical factors that may influence use of serum prostate specific antigen and its isoforms in a screening programme for prostate cancer. Ann Clin Biochem (38) 188-216
- Ferlay J et al. (2013) Cancer incidence and mortality patterns in Europe: Estimates for 40 countries in 2012. Europ J Cancer (49) 1374-1403

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14. Referências

1. Milford Ward A. et al. (2001) Prostate specific antigen: biology, biochemistry and available commercial assays. Ann Clin Biochem 38: 633-651
2. Ballentine Carter H. (2004) Prostate Cancers in Men with Low PSA Levels – Must We Find Them? N Engl J Med 350(22): 2292-2294
3. Jain S et al. (2002) Improving the utility of prostate specific antigen (PSA) in the diagnosis of prostate cancer: the use of PSA derivatives and novel markers. Postgrad Med J (78) 646-650
4. Price et al. (2001) Pre- and post-analytical factors that may influence use of serum prostate specific antigen and its isoforms in a screening programme for prostate cancer. Ann Clin Biochem (38) 188-216
5. Ferlay J et al. (2013) Cancer incidence and mortality patterns in Europe: Estimates for 40 countries in 2012. Europ J Cancer (49) 1374-1403
6. Leitlinienprogramm Onkologie (Deutsche Krebsgesellschaft, Deutsche Krebshilfe, AWMF): Interdisziplinäre Leitlinie der Qualität S3 zur Früherkennung, Diagnose und

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nal von minden GmbH